CRP Testkit

Deutsch

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English

For human medical use only!

for quantitative in vitro determination of c-reactive protein in human finger blood, EDTA venous blood or serum on the smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0100	CRP test kit	32 tests
Order number: ST1000	CRP control kit	2 x 2 ml (low/high)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. It increases already after 6 hours in acute inflammatory processes associated with bacterial infections, post operative conditions or tissue damage. Many studies have demonstrated that in apparently healthy subjects there is a direct correlation between CRP concentrations and the risk for developing coronary heart disease (CHD).

Method

Measurement at 546 nm or 700 nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

Measurement Range

Immunturbidimetric Assay:

Serum: 0.5 - 120 mg/l at 546 nm / 1.0 - 120 mg/l at 700 nm

Whole blood: 2.0 - 240 mg/l at 546 nm & 700 nm (LOT and HCT pending)

Stability and Storage

Storage at 2 - 8 °C

Unopened vacuum package: Until expiration date

Opened vacuum package: 3 months

Storage at 18 - 25 °C

Unrefrigerated: 1 month

Refrigerated at least 16h/day: 2 months

Warnings and Precautions

The reagents contain sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Freshly drawn fingertip blood is preferred. Alternatively EDTA venous blood or serum can be used. *EDTA plasma samples show ~16% lower values, Li-Hep plasma samples show ~28% lower values compared to serum samples.* After pipetting the sample into the ERS cuvette it must be processed with a CUBE/smart laboratory photometer within 4 hours.

Reference Range

Values < 10 mg/l are expected as normal. Any condition that results in sudden or severe inflammation may increase CRP level. It is recommended that each laboratory establish its own expected range.

Quality Control

For internal quality control the CRP control kit should be used.

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