# Cystatin C (GFR) test kit

for quantitative in vitro determination of Cystatin C in human capillary blood or serum on the smart or CUBE laboratory photometer.

Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

For human medical use only!

### Order information

Order number: ST0240  
Order number: ST2400

### Indication

Cystatin C (GFR) test kit  
Cystatin C (GFR) control kit

### Kit size

32 tests  
2 x 1 ml (low/high)

**Test kit preparation:** Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

**Reagent appearance in ERS cuvette has to be clear and colourless. If a colour change to yellowish or brownish colour occurs, the respective ERS cuvette should not be used!**

## Summary

Cystatin C is a proteinase inhibitor with a low molecular mass of 13Kda that is produced at a constant rate in all nucleated cells and appears in human plasma and serum. Cystatin C is freely filtered through the glomerulus, is not secreted by the tubule or eliminated via any extra-renal route, and is almost completely absorbed and catabolized by proximal tubular cells. Therefore, the plasma concentration of Cystatin C is almost exclusively determined by the glomerular filtration rate (GFR), making Cystatin C an excellent indicator of GFR. A number of clinical studies have shown that Cystatin C is more accurate than plasma creatinine and the Cockcroft-Gault estimation of creatinine clearance and is more reliable than the 24-h creatinine clearance. There is a growing body of evidence that suggests that Cystatin C can be used to detect kidney disease at earlier stages than serum creatinine which may help facilitate prevention efforts in the elderly and those with diabetes, hypertension, or cardiovascular disease.

## Method / Measurement Range

Cystatin C concentration: 0.33 - 7.0 mg/l  
Equals GFR (glomerular filtration rate):  
541 - 3 (ml/min/1.73 m<sup>2</sup>)

Note: The measurement range is LOT pending!

## Principle

Turbidimetric measurement of the Cystatin C concentration at 700 nm wave length.

## Test Kit

ERS cuvette filled with: Tris-buffer solution  
ERS cap filled with: Cystatin C antibody reagent

## Stability and Storage

Stable when stored at 2 - 8 °C until the expiration date stated on the labels. DO NOT FREEZE.

## Warnings and Precautions

This test kit is for in vitro diagnostic use only. Specimens containing human sourced materials should be handled as if potentially infectious using safe laboratory procedures. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Avoid ingestion and contact with skin and eyes.

## Waste Management

Please refer to local legal requirements.

## Specimen Sample Material

20 µl fresh capillary blood from fingertip or serum.

## Reference Range

See „Reference Graph“. However, it is recommended to check the range for the regional population.

## Quality Control

For internal QC use the Cystatin C (GFR) control kit only.  
Order no.: ST2400

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ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 -25 °C) before use!

Correct use of the capillary

- A Insert the plunger into the capillary opening at the **black marker**. Then, hold capillary in a sloping position to blood drop until it is completely filled up to the **red marker**.  
**Aspirate**
- B WRONG: Capillary is not completely filled with blood up to the **red marker**. Hold the capillary horizontally when aspirating to ensure it gets filled completely.
- C CORRECT: Capillary is completeley filled with blood up to the **red marker**  
**Dispense**
- D WRONG: Plunger has not been pushed all the way into the capillary. Blood is left in capillary.
- E CORRECT: Plunger has been pushed completely into the capillary. **White stopper** is pressed to the end of the capillary. No blood is left.

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touchscreen

2. Select whole blood or serum

- 2.1 Sampletype menu: Select the sample type
- 2.2a Haematocrit correction (if desired)
- 2.2b Automatic haematocrit correction (smart 700/546 and CUBE-S only!)
- 2.3a Aspirate 20 µl fresh capillary blood from fingertip using the capillary  
OR ...
- 2.3b Aspirate 20 µl serum from primary tube using the pipette
- 2.4a/b Dispense 20 µl sample into ERS cuvette INTO LIQUID.

3. Sample processing

- 3.1 Apply ERS cap firmly onto ERS cuvette
- 3.2 Place ERS cartridge into laboratory photometer
- 3.3 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.