# D-Dimer Control Kit

## Deutsch

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## English

For human medical use only!

for quality control procedures of the quantitative in vitro determination with the D-Dimer test kit on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0150	D-Dimer test kit	16 tests
Order number: ST0151	D-Dimer test kit	6 tests
Order number: ST1500	D-Dimer control kit	2 x 1 ml (low/high)

**Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

### Summary

This D-dimer control kit is intended for use in the clinical laboratory as a control for quantitative determination of D-dimer in the normal range and high range. Controls are prepared from human blood and each single plasma used for this lots are tested HBsAg, HIV 1/2 Ab and HCV Ab negative. However, the samples should be handled according to general rules in a clinical laboratory environment. This control kit is for in vitro diagnostic (IVD) use only.

### Storage and Handling

Store at 2 - 8 °C. Kit can be used until expiry date printed on the labels. DO NOT FREEZE!

### After reconstitution:

Stored at	Expiration
2 - 8 °C	14 days

### Waste Management

Please refer to local legal requirements. All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

### Preparation

The D-dimer controls are lyophilized and have to be reconstituted with 1 ml of distilled water 30 minutes prior to use. Note: Once reconstituted there is no further preparation step necessary. (such as spinning or mixing with citrate) It can directly be used same way like the patient sample!

### Control Values and Calibrator Standardization

No international calibrator is available for the standardization of D-dimer. Therefore the D-dimer calibrator is calibrated by using a commercial ELISA kit and the control values are assigned against the master calibration ELISA-based. Accuracy should be checked according to your local/national QC requirements with D-dimer control high and D-dimer control low. The control plasma is treated in the same way as patient samples and results should be checked against the acceptable confidence limits supplied on the kit insert sheet and the vial label.

### Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Measure controls the same way like a patient sample and compare the values indicated in the lot-specific value sheet included in the control kit.

If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

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