# Lipoprotein (a) test kit

for quantitative in vitro determination of lipoprotein (a) on smart or CUBE laboratory photometer.

Eurolyser Diagnostica GmbH  
Bindergrasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order Information

Order number: ST 0141  
Order number: ST 1400

### Indication

Lipoprotein (a) test kit  
Lipoprotein (a) control kit

### Kit size

16 tests/kit  
2 x 1 ml (low/high)

**Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

Lipoprotein (a), Lp(a), is a low density lipoprotein like particle containing apolipoprotein B-100 disulphide-linked to the large glycoprotein called apolipoprotein (a). The characteristic feature of Lp(a) is that it is distinct from all other serum proteins and apolipoproteins. This protein is believed to be inherited as an autosomal dominant trait and appears to be insensitive to either diet, lifestyle or most hypolipidaemic drugs. Since its discovery by BERG in 1963 there has been a considerable rise in interest not only in specialised research centers but also in clinical laboratories in the accurate measurement of Lp(a) in blood.

This interest was stimulated by reports indicating that levels above 20 - 30 mg/dl present in approx. 25% of the population are associated with an increased risk of coronary heart disease. Many investigators have confirmed that high Lp(a) concentrations represent an indicator of risk for cardiovascular disease, especially when serum LDL or APO B are elevated.

Therefore a convenient and reliable method for the quantification of Lp(a) in serum or plasma is important for identification of individuals at risk for developing arteriosclerosis.

## Method

Immunoturbidimetric kinetic measurement of the Lp(a) concentration at 700 nm (Absorbance).

## Measurement Range

Measurement range: 0 - 100 mg/dl (LOT-pending)

Samples with concentrations higher than the upper limit of the measurement range must be diluted 1 + 2 with physiological saline (0.9% NaCl solution), e.g. 10 µl sample + 20 µl 0.9% NaCl solution, and the result multiplied by 3.

## Sample Material

Fresh sample material: Serum or EDTA plasma. Immediately centrifuge and process after blood collecting. Very lipaemic or turbid specimen must be clarified before assay is performed.

## Reference Range

for adults < 30 mg/dl.

Each laboratory is recommended to establish a range of normal values for the population in their region.

## Test Kit

ERS cuvette filled with glycine buffer solution

ERS cap filled with latex reagent - anti Lp(a) antibodies

## Quality Control

For internal quality control the Lipoprotein (a) control kit is recommended. Order number: ST 1400

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

## Warning and Precautions

This test kit is for in vitro diagnostic use only!

DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Observe all necessary precautions for the use of laboratory reagents.

## Waste Management

Always refer to local legal requirements.

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## ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press "Measurement" button, enter required information using the touchscreen

### 2. Sample preparation

For further details please refer to the user manual of the laboratory photometer.

- 2.1 Aspirate 10 µl sample from centrifuged sample tube
- 2.2 Dispense 10 µl sample into ERS cuvette INTO liquid

### 3. Test preparation

- 3.1 Apply ERS cap
- 3.2 Place ERS cartridge into laboratory photometer
- 3.3 Start automatic sample processing by pressing the  start button on the smart laboratory photometer or by closing the door of the CUBE laboratory photometer.