# NH3 Ammonia test kit

for quantitative in vitro determination of ammonia on smart laboratory photometer.

Eurolyser Diagnostics GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order Information

Order number: ST0500  
Order number: ST5000

IVD

### Indication

NH3 Ammonia test kit  
NH3 Ammonia control kit

### Kit size

16 tests/kit  
1 x 5ml (dec level)

**Test kit preparation: Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

The bulk of ammonia in the body is generated in the gastrointestinal system by action of bacterial enzymes on the contents of the colon and from hydrolysis of glutamine. It is removed in the liver and converted to urea through a series of enzymatic reactions in the Krebs-Henseleit cycle. Among other conditions, advanced liver disease and hepatic encephalopathy result in elevated levels of ammonia in blood. Hyperammonemia is also common in inherited deficiencies of the enzymes involved in the conversion of ammonia to urea. The determination of ammonia is very useful in the diagnosis and prognosis of Reye's Syndrome. Elevated blood ammonia exerts toxic effects on the central nervous system.

## Principle

Ammonia reacts with alpha-ketoglutarate and reduced cofactor to form L-glutamate and the cofactor. The reaction is catalyzed by glutamate dehydrogenase. The decrease in absorbance due to the oxidation of the reduced cofactor is monitored at 340nm and is proportional to the ammonia concentration.

## Sample Material

The accuracy of ammonia determination is extremely dependent on sample collection. Use only fresh EDTA plasma (Ammonia Heparin Plasma or Serum must not be used). After drawing the blood, immediately centrifuge the collection tube to get plasma. The plasma sample has to be analysed within 5 minutes on the analyser.

## Typical Measurement Range

Measurement range\*: 10 - 1100 µmol/l  
Lower detection limit\*: 10 µmol/l

\* Lot depending

## Reference Range

Adults: 12 - 47 µmol/l

Each laboratory, however, is recommended to establish a range of normal values for the population in their region.

## Testkit

ERS Cuvette pre-filled with: NH3 reagent

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

## Quality Control

For internal quality control the NH3 Ammonia control kit should be used. Order number: ST5000.

## Warnings and Precautions

This testkit is for in vitro diagnostic use only. Contains sodium azide (0,95 g/l). Do NOT SWALLOW! Avoid contact with skin and eyes. Take the necessary precautions for the use of laboratory reagents.

## Waste Management

Please refer to local legal requirements.

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---

## ATTENTION!

Allow single test min. 10 minutes to warm to room temperature!

### 1. Preparation of test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Pull of ERS cuvettecap
- 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser

---

### 2. Sample preparation plasma

- 2.1 Aspirate 20 µl plasma from primary tube

---

### 3. Test preparation

- 3.1 Pipette 20 µl sample into ERS cuvette
- 3.2 Apply ERS cuvettecap
- 3.3 Place ERS cuvette into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the smart laboratory photometer.