# NH<sub>3</sub> Ammonia test kit

for quantitative in vitro determination of ammonia on a smart photometer.

Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order information

Order number: ST0500  
Order number: ST5000

### Indication

NH<sub>3</sub> Ammonia test kit  
NH<sub>3</sub> Ammonia control kit

### Kit size

16 tests/kit  
1 x 5 ml (decision level)

 **Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

The bulk of ammonia in the body is generated in the gastrointestinal system by action of bacterial enzymes on the contents of the colon and from hydrolysis of glutamine. It is removed in the liver and converted to urea through a series of enzymatic reactions in the Krebs-Henseleit cycle. Among other conditions, advanced liver disease and hepatic encephalopathy result in elevated levels of ammonia in blood. Hyperammonaemia is also common in inherited deficiencies of the enzymes involved in the conversion of ammonia to urea. The determination of ammonia is very useful in the diagnosis and prognosis of Reye's Syndrome. Elevated blood ammonia exerts toxic effects on the central nervous system.

## Method

Ammonia reacts with alpha-ketoglutarate and reduced co-factor to form L-glutamate and the cofactor. The reaction is catalyzed by glutamate dehydrogenase. The decrease in absorbance due to the oxidation of the reduced cofactor is monitored at 340 nm and is proportional to the ammonia concentration.

## Test Kit

ERS Cuvette pre-filled with: NH<sub>3</sub> reagent

## Typical Measurement Range

Measurement range\*: 10 - 1100 µmol/l  
Lower detection limit\*: 10 µmol/l

\*LOT specific

## Sample Material

40 µl of fresh EDTA plasma is needed. The accuracy of ammonia determination is extremely dependent on sample collection. Ammonia Heparin Plasma or Serum must not be used. After drawing the blood, immediately centrifuge the collection tube to get plasma. The plasma sample has to be analysed within 5 minutes on the analyser.

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

## Reference Range

Adults: 12 - 47 µmol/l

It is recommended that every laboratory establishes its own reference ranges.

## Quality Control

For internal quality control the Eurolyser NH<sub>3</sub> Ammonia control kit should be used. Order number: ST5000

## Precision

N = 20; CV (246 µmol/l) < 4.9%

## Correlation

N = 20; R<sup>2</sup> = 0.985; y = 0.972x + 0.233  
y = Eurolyser NH<sub>3</sub>, x = reference method NH<sub>3</sub>

## Warnings and Precautions

This test kit is for in vitro diagnostic use only. Contains sodium azide (0,95 g/l). **DO NOT SWALLOW!** Avoid contact with skin and eyes. Take the necessary precautions for the use of laboratory reagents.

## Waste Management

Please refer to local legal requirements.

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