ASO (anti-Streptolysin-O) Testkit

Deutsch

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Literatur / References:

1. Klein G.C., Baker C.N., Jones W.L., Upper Limits of normal Antistreptolysin-O and Antideoxyribonuclease B titres. Applied Microbiology. 1971 ; 21: 999-1001
2. Thomas L. ed., Clinical Laboratory Diagnosis.Streptococcus pyrogenes Infection. In: Thomas L.ed., Clinical Laboratory Diagnosis. Use and assessment of Clinical Laboratory results. 1st Edition (1998). TH-Books, Frankfurt/Main Germany.

ASO (anti-Streptolysin-O) test kit

English

for quantitative in vitro determination of anti-streptolysin-O on the smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0250	ASO test kit	16 tests/kit
Order number: ST2500	ASO control kit	1 x 1 ml (decision range)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The ASO test is used to determine post streptococcal complications including rheumatic fever and glomerulonephritis as well as recent streptococcal infection. The presence and level of ASO antibodies in human blood/serum directly reflects the extent and degree of infection. Elevated levels of ASO may also be present in other conditions including tonsillitis, acute rheumatoid arthritis, scarlet fever, and various other streptococcal infections.

Method / Measurement Range

Latex-enhanced immunturbidimetric assay: 100 - 2000 IU/ml

Principle

Latex particles coated with Streptolysin-O react with ASO in the patient sample resulting in an increase in turbidity. The absorbance is measured at 700 nm and is directly proportional to the concentration of ASO.

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

The reagents contain sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Freshly drawn whole blood from fingertip is preferred. Alternatively EDTA venous blood or serum can be used.

Reference Range

≤ 100 IU/ml for preschool children
≤ 250 IU/ml for school-age children
≤ 200 IU/ml for adults
It is recommended that each laboratory should establish its own expected range.

Quality Control

For internal quality control the ASO control kit should be used.
Order number: ST2500

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ATTENTION!

Use the integrated capillary ONLY for aspirating whole blood from finger tip! Add EDTA venous blood or serum ONLY with pipette! Allow single test at least 10 minutes to warm up to room temperature before use!

Correct use of the integrated capillary

- A CORRECT: Capillary is completely filled with blood
- B WRONG: Capillary is partially filled with blood
- C WRONG: Too much blood around the capillary and/or ERS cap

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- 1. Preparation of test system**
- 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Place ERS cap in test kit rack
 - 1.4 Press „Measurement“ button, enter required information using the touch screen

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
- 2. Selection of whole blood or serum**
- Please consult the laboratory photometer user manual for more information.
- 2.1 Sampletype menu: Select the sample type
 - 2.2a Haematocrit correction (if desired)
 - 2.2b Automatic haematocrit-correction (smart 700/546 and CUBE-S only!)
 - 2.3 Instrument ready for analysis

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- 3. Sample preparation with whole blood from fingertip**
- 3.1a Aspirate 5 µl whole blood with integrated capillary from finger tip
 - 3.2a Apply ERS cap firmly onto ERS cuvette
 - 3.3a Shake ERS cartridge 3 - 4x gently until the finger blood is fully dispensed out of capillary into ERS liquid
 - OR with EDTA venous blood or serum**
 - 3.1b Aspirate 5 µl EDTA venous blood or serum with pipette
 - 3.2b Dispense 5 µl EDTA venous blood or serum INTO LIQUID in ERS cuvette
 - 3.3b Apply ERS cap firmly onto ERS cuvette
 - 3.4 Place ERS cartridge into laboratory photometer
 - 3.5 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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Technical details subject to change without notice.