# CRP Testkit

## Deutsch

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# CRP test kit

## English

for quantitative in vitro determination of c-reactive protein in finger blood, EDTA venous blood or serum on the smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0100	CRP test kit	32 tests/pack
Order number: ST0102	CRP test kit with integr. capillary	32 tests/pack
Order number: ST1000	CRP control kit	2 x 2 ml (low/high)

**Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

### Summary

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. It increases already after 6 hours in acute inflammatory processes associated with bacterial infections, post operative conditions or tissue damage. Many studies have demonstrated that in apparently healthy subjects there is a direct correlation between CRP concentrations and the risk for developing coronary heart disease (CHD).

### Method

Measurement at 546 nm or 700 nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

### Measurement Range

Immunturbidimetric Assay:  
 Serum: 0.5 - 120 mg/l at 546 nm / 1.0 - 120 mg/l at 700 nm  
 Whole blood: 2.0 - 240 mg/l at 546 nm & 700 nm (LOT and HCT pending)

### Stability and Storage

**Storage at 2 - 8 °C**  
 Unopened vacuum package: Until expiration date  
 Opened vacuum package: 3 months  
**Storage at 18 - 25 °C**  
 Unrefrigerated: 1 month  
 Refrigerated at least 16h/day: 2 months

### Warnings and Precautions

The reagents contain sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

### Waste Management

Please refer to local legal requirements.

### Specimen Sample Material

Freshly drawn fingertip blood is preferred. Alternatively EDTA venous blood or serum can be used. *EDTA plasma samples show ~16% lower values, Li-Hep plasma samples show ~28% lower values compared to serum samples.* After pipetting the sample into the ERS cuvette it must be processed with a CUBE/smart laboratory photometer within 4 hours.

### Reference Range

Values < 10 mg/l are expected as normal. Any condition that results in sudden or severe inflammation may increase CRP level. It is recommended that each laboratory establish its own expected range.

### Quality Control

For internal quality control the CRP control kit should be used.  
 Order number: ST1000

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