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CRP test kit with integrated capillary

English

for quantitative in vitro determination of c-reactive protein in finger blood, EDTA venous blood or serum on the smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0100	CRP test kit	32 tests/pack
Order number: ST0102	CRP test kit with integr. capillary	32 tests/pack
Order number: ST1000	CRP control kit	2 x 2 ml (low/high)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. It increases already after 6 hours in acute inflammatory processes associated with bacterial infections, post operative conditions or tissue damage. Many studies have demonstrated that in apparently healthy subjects there is a direct correlation between CRP concentrations and the risk for developing coronary heart disease (CHD).

Method

Measurement at 546 nm or 700 nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

Measurement Range

Immunturbidimetric Assay:
 Serum: 0.5 - 120 mg/l at 546 nm / 1.0 - 120 mg/l at 700 nm
 Whole blood: 2.0 - 240 mg/l at 546 nm & 700 nm (LOT and HCT pending)

Stability and Storage

Storage at 2 - 8 °C
 Unopened vacuum package: Until expiration date
 Opened vacuum package: 3 months
Storage at 18 - 25 °C
 Unrefrigerated: 1 month
 Refrigerated at least 16h/day: 2 months

Warnings and Precautions

The reagents contain sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Freshly drawn fingertip blood is preferred. Alternatively EDTA venous blood or serum can be used. *EDTA plasma samples show ~16% lower values, Li-Hep plasma samples show ~28% lower values compared to serum samples.* After pipetting the sample into the ERS cuvette it must be processed with a CUBE/smart laboratory photometer within 4 hours.

Reference Range

Values < 10 mg/l are expected as normal. Any condition that results in sudden or severe inflammation may increase CRP level. It is recommended that each laboratory establish its own expected range.

Quality Control

For internal quality control the CRP control kit should be used.
 Order number: ST1000

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ATTENTION!

Use the integrated capillary ONLY for aspirating whole blood from finger tip! Add EDTA venous blood or serum ONLY with pipette! Allow single test at least 10 minutes to warm up to room temperature before use!

Correct use of the integrated capillary

- A CORRECT: Capillary is completely filled with blood
- B WRONG: Capillary is partially filled with blood
- C WRONG: Too much blood around the capillary and/or ERS cap

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2. Selection of whole blood or serum

Please consult the laboratory photometer user manual for more information.

- 2.1 Samplotype menu: Select the sample type
- 2.2a Haematocrit correction (if desired)
- 2.2b Automatic haematocrit-correction (smart 700/546 and CUBE-S only!)
- 2.3 Instrument ready for analysis

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3. Sample preparation with whole blood from fingertip

- 3.1a Aspirate 5 µl whole blood with integrated capillary from finger tip
- 3.2a Apply ERS cap firmly onto ERS cuvette
- 3.3a Shake ERS cartridge 3 - 4x gently until the finger blood is fully dispensed out of capillary into ERS liquid

OR with EDTA venous blood or serum

- 3.1b Aspirate 5 µl EDTA venous blood or serum with pipette
- 3.2b Dispense 5 µl EDTA venous blood or serum INTO LIQUID in ERS cuvette
- 3.3b Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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
Technical details subject to change without notice.