D-Dimer Testkit

Deutsch

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D-Dimer test kit

English

for quantitative in vitro determination of D-dimer on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0150	D-Dimer test kit	16 tests/pack
Order number: ST0151	D-Dimer test kit	6 tests/pack
Order number: ST1500	D-Dimer control kit	2 x 1 ml (high/low)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Elevated levels of D-dimer are found in patients with confirmed deep venous thrombosis (DVT), pulmonary embolism (PE), DIC, and trauma. D-dimer levels rise during pregnancy and high levels are associated with complications.

Limitation and Interferences

Each measurement series should be checked with D-dimer control high and low. Presence of rheumatoid arthritis factor may result in false-positive results (influence not quantified). No interference is found from CRP (< 50 mg/l), bilirubin (< 0.27 g/l), haemoglobin (< 6.7 g/l), or triglycerides (< 3.3 g/l). Results from patients with heterophilic antibody should be interpreted with caution since this test kit contains mouse antibodies and interference may occur resulting in falsely elevated or decreased values.

Method

Kinetic measurement of the D-dimer concentration at 700 nm (absorbance)

Measurement Range

500 - 8500 ng/ml (LOT pending)

Test Kit

ERS cuvette filled with buffer
ERS cap filled with antibodies (Mouse JIF-23) latex reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Sample Material

Citrat plasma (9+1 ratio) is the first choice (Set samplotype: „Cit. Plasma“) Alternatively, EDTA-plasma can be used (Set samplotype: „Plasma“)

Reference Range

Samples with a D-dimer concentration > 500 ng/ml are considered as pathologic.

Quality Control

For internal quality control the D-Dimer control kit should be used.
Order number: ST1500

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