# iFOB Testkit

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## English

For human medical use only!

for quantitative-immunological in vitro determination of faecal occult blood on the smart or CUBE laboratory photometer.

## Order information

Order number: ST0200  
Order number: ST0205  
Order number: ST2000

## Indication

iFOB test kit  
iFOB sample collector kit  
iFOB control kit

## Kit size

32 tests  
32 collectors  
2 x 2 ml (low/high)

**Test kit preparation:** Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

## Summary

Haemoglobin is the major protein in red blood cells. The presence of haemoglobin in feces is indicative for the presence of blood (Faecal Occult Blood, FOB) as a result of bleedings associated to pathologies of the gastrointestinal tract such as colon polyps, adenomas, colorectal carcinoma, ulcerative colitis and Crohn's disease. Colorectal cancer (CRC) is the third most common form of cancer in the western world, however, it often causes no symptoms until it has reached a relatively advanced stage. Thus, many organisations recommend the screening of this disease through iFOB testing. Colorectal bleeding can be sometimes intermittent, therefore negative results do not necessarily rule out the disease; a colonoscopy is recommended as a definitive test! iFOB testing is also not suitable as screening method for diseases of the upper gastrointestinal tract.

## Method

Polyclonal antibody Hb detection in faeces. The latex particles coated with antibodies anti-human haemoglobin agglutinate when they react with faeces samples containing human haemoglobin. Photometric measurement at 700 nm of the turbidimetric complex.

## Measurement Range

25 - 1000 µg/l (LOT dependent)

Samples exceeding 1000 µg/l should be diluted 1 + 5 with physiological saline (0.9%-NaCl), the result multiplied by 6 and then retested.

## Sample Material

For collecting a faeces sample use the iFOB sample collector kit exclusively and follow the instructions provided with the kit. Order number: ST0205

## Test kit

ERS cap filled with latex particles coated with polyclonal anti human haemoglobin antibodies, pH 8.2.

ERS cuvette filled with glycine buffer 100 mmol/l, pH 10

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2 - 8 °C) from the date of opening. DO NOT FREEZE!

## Warning and Precautions

DO NOT INGEST! This test kit is for in vitro diagnostic use only. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds.

## Waste Management

Please refer to local legal requirements.

## Reference Range

It is recommended that each laboratory establishes its own cut-off according to the local population and the purpose of the assay (e.g. screening, follow-up, etc.).

## Sensitivity and specificity vary depending on the cut-off as follows:

Cut-off (µg/l)	Sensitivity (%)	Specificity (%)
50	79,4	89,7
100	76,5	95,3
150	70,6	95,9

## Quality Control

For internal quality control the iFOB control kit is recommended.  
Order number: ST2000

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