# Haemoglobin Control Kit

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# Haemoglobin Control Kit

for quality control procedures of the quantitative in vitro determination with the haemoglobin test kit on smart or CUBE (-S) laboratory photometer.

Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order information

Order number: ST0190  
Order number: ST0195  
Order number: ST1900

IVD

### Indication

Haemoglobin test kit  
Haemoglobin test kit with capillary  
Haemoglobin control kit

### Kit size

32 Tests/kit  
32 Tests/kit  
2 x 1.5 ml (low/high)

**Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

The haemoglobin control kit is intended for use in the clinical laboratory as a control for quantitative determination of HGB in the normal range and high range. It is prepared from haemoglobin equivalent.

## Preparation and Handling

The control is liquid and ready to use.

Thoroughly mix the content of the vial before each use by gently inverting it. Turbidity shows bacterial decontamination, the vial has to be discarded. Always use disposable pipette tips!

## Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date on the label. Once opened vials are stable until the expiry date when stored tightly capped at 2 - 8 °C.

## Waste Management

Please refer to local legal requirements.

## Warnings and Precautions

For in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Set the sex to „Control“ before running the QC.

Compare the results with the values indicated in the lot-specific value table.


If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

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
## Processing of a Haemoglobin control measurement

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


## English

### ATTENTION!

Allow single test to warm up at room temperature for a minimum of 10 minutes before use!

Test can only be processed on instruments with firmware version **1.33** or higher!

#### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in the test kit rack
- 1.3 Detach ERS cap
- 1.4 Press „Measurement“ button, enter required information using the touch screen

#### 2. Photometer Settings

- 2.1 Sex menu: Select „Control“
- 2.2 Sampletype menu: Select „Pipette 20µ“
- 2.3 Instrument ready for analysis

#### 3. Processing control measurement with kit ST0195

- 3.1 Dispense some drops control liquid into serum cup and aspirate 20 µl using a pipette.
- 3.2 Dispense 20 µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer and start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE(-S) laboratory photometer.

#### 4. Processing control measurement with kit ST0190

- 4.1 Dispense some drops control liquid into serum cup and aspirate 20 µl from serum cup.
- 4.2 Dispense 20 µl sample INTO LIQUID of ERS cuvette
- 4.3 Apply ERS cap firmly onto ERS cuvette
- 4.4 Gently shake the ERS cartridge 3 to 4 times.
- 4.5 Place ERS cartridge into laboratory photometer and start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

## Wertetabelle HGB control kit

## - Value sheet HGB control kit

LOT

### Niedrige Kontrolle - Low control

Ziel / Target	Wert / Value (min.)	Wert / Value (max)

### Hohe Kontrolle - High control

Ziel / Target	Wert / Value (min.)	Wert / Value (max)

Technical details subject to change without notice.

Doc-ID: e2689b5c-da26-4709-a944-89f230e92593 ; latest revision: DP6.0 ; date of print: 2019-03-11