# smart HGB test kit

Testkit for quantitative determination of Hemoglobin from whole blood on smart photometer.

M Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order Information

Order No.: ST0190  
Order No.: ST1900

### Indication

smart HGB test kit  
Haemoglobin control kit

### Kit size

32 tests/kit  
2 x 1.5 ml (low / high)

## Summary

In routine medical diagnostics, the test may be used for the following purposes: To diagnose anemia or polyglobulia, to monitor risk groups for iron deficiency (pregnant women, infants, blood donors, hemodialysis patients, diabetics, athletes, ...), as part of point of care testing when acute bleedings or hemorrhages are suspected, at parasitic diseases as well as quantifying blood loss.

## Assay principle

The reagent used binds hemoglobin and its variants to form a coloured complex, which is detected photometrically at 546nm.

Measuring Range: 0-30 g/dl

Lower Detection Limit: 0.5 g/dl

## Testkit

Filled with AHD reagent.

## Stability and Storage

Store at 2-40°C. DO NOT FREEZE. Kit can be used until expiry date printed on the labels.

## Warning and Precautions

This testkit is for in vitro diagnostic use only. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds. Take the necessary precautions for the use of laboratory reagents. Do NOT INGEST!

## Waste Management

Please refer to local legal requirements.

## Sample Material

Fresh capillary blood is preferred.  
Alternatively whole blood from EDTA can be used\*.  
Make sure sample is mixed well.

\*The assay is calibrated for capillary finger blood, please verify your results if venous whole blood is used and set a correlation factor/offset according to your validation.

## Reference Ranges

Women	12-16 g/dl
Men	14-18 g/dl
Infants	10-15 g/dl
Children	12-16 g/dl

This range is given for orientation only. Each Laboratory should establish his own reference values.

## Control Material

For internal quality assurance the Haemoglobin control kit should be used; Order No.: ST1900

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## ATTENTION!

Allow single test min. 10 minutes to warm to room temperature (20-25 °C)!  
Test can be processed on instruments with firmware version 1.33 or higher only.

## English

### 1. Preparation of test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette in smart test kit rack
- 1.3 Pull off ERS cuvettecap
- 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser

### 2. Sample preparation for whole blood (fingertip or EDTA primary tube)

- 2.1a Aspirate 20µl fresh capillary blood from fingertip OR ...
- 2.1b Aspirate 20µl whole blood from EDTA primary tube
- 2.2 Pipette 20µl sample INTO LIQUID in ERS cuvette
- 2.3 Apply ERS cuvettecap firmly onto ERS cuvette again

### 3. Test preparation

- 3.1 Shake closed ERS cuvette 6-8x vigorously
- 3.2 Place ERS cuvette into laboratory photometer
- 3.3 Start automatic sample processing by pressing the  start button