Haemoglobin test kit

with capillary

for quantitative in vitro determination of haemoglobin from whole blood on smart or CUBE-S laboratory photometer.

EUROLyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order Information

Order No.: ST0195
Order No.: ST1900

Indication

Haemoglobin test kit with capillary
Haemoglobin control kit

Kit size

32 tests/kit
2 x 1,5 ml (low / high)

Summary

In routine medical diagnostics, the test may be used for the following purposes: To diagnose anaemia or polycythemia, to monitor risk groups for iron deficiency (pregnant women, infants, blood donors, haemodialysis patients, diabetics, athletes, ...), as part of point-of-care testing when acute bleedings or haemorrhages are suspected, for detecting parasitic diseases as well as for quantifying blood loss.

Method / Measurement Range

The reagent used binds haemoglobin and its variants to form a coloured complex, which is detected photometrically at 546nm.

Measurement Range: 0 - 30 g/dl

Lower Detection Limit: 0.5 g/dl

Testkit

Filled with AHD reagent.

Stability and Storage

Store at 2 - 40°C. Kit can be used until expiry date printed on the labels. DO NOT FREEZE!

Warnings and Precautions

This testkit is for in vitro diagnostic use only. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents. DO NOT INGEST!

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Fresh capillary blood is preferred.

Alternatively whole blood from EDTA can be used*.

Make sure sample is mixed well.

*The assay is calibrated for capillary finger blood. Please verify your results if venous whole blood is used and set a correlation factor/offset according to your validation.

Reference Ranges

Women : 12 - 16 g/dl

Men: 14 - 18 g/dl

Infants: 10 - 15 g/dl

Children: 12 - 16 g/dl

This range is given for orientation only. Each laboratory should establish its own reference values.

Quality Control

For internal quality assurance the Haemoglobin control kit should be used; Order No.: ST1900

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ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C)!

Test can be processed on instruments with firmware version 1.33 or higher only.

Correct use of capillary

- A Insert the plunger into the capillary opening at the **wide, black marker**. Then hold capillary in a sloping position to blood drop until it is filled up completely to the small red marker.
Aspirate
- B INCORRECT: Capillary is not filled up completely to the small red marker with blood. Keep capillary horizontal when aspirating to ensure it gets filled completely.
- C CORRECT: Capillary is completely filled up with blood to the small red marker.
Dispense
- D INCORRECT: Plunger has not been completely pressed into capillary. Blood is left in capillary.
- E CORRECT: Plunger has been completely pressed into capillary. White stopper is pressed to the end of the capillary. No blood is left.

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Pull off ERS cuvette cap
- 1.4 Press symbol on touch screen, required information is entered over the analyser's touch display

2. Sample preparation for whole blood (fingertip or EDTA primary tube)

- 2.1a Choose sample type 'CapilBlood' if sample is aspirated from fingertip
- 2.2a Aspirate 20µl from fingertip (See "Correct use of capillary"!)
OR ...
- 2.1b Choose sample type 'Pipette20µl' if sample is aspirated from EDTA primary tube
- 2.2b Aspirate 20µl whole blood from EDTA primary tube

3. Test preparation

- 3.1 Pipette 20µl sample INTO LIQUID in ERS cuvette
- 3.2 Apply ERS cuvette cap firmly onto ERS cuvette again
- 3.3 Place ERS cuvette into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the smart photometer or by closing the door of the CUBE-S photometer