Homocysteine test kit

for quantitative in vitro determination of homocysteine in serum or plasma on the smart laboratory photometer.

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

For human medical use only!

Order information

Order number: ST0120
Order number: ST1200
Order number: ST1210

Indication

Homocysteine test kit
Homocysteine control kit
Homocysteine calibrator kit

Kit size

16 tests
2 x 3 ml (low/high)
1 x 1 ml

 **Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

Summary

Homocysteine is a thiol containing amino acid produced by the intracellular demethylation of methionine. Total homocysteine (tHCY) represents the sum of all forms of homocysteine including forms of oxidized, protein-bound and free forms. Elevated level of tHCY has emerged as an important risk factor in the assessment of cardiovascular disease. Excess homocysteine in the bloodstream may cause injuries to arterial vessels due to its irritant nature and result in inflammation and plaque formation, which may eventually cause blockage of bloodflow to the heart. Elevated tHCY levels are caused by four major factors:

- Genetic deficiencies
- Nutritional deficiency in B vitamins such as B6, B12 and folate
- Renal failure for effective amino acid clearance
- Drug interactions such as nitricide, methotrexate and phenytoin

Elevated levels of tHCY are also linked with Alzheimer's disease and osteoporosis

Method

Kinetic determination based on an enzymatic cycling method

Measurement Range

Serum/Plasma: 3 - 50 µmol/l (lot specific)

Test Kit

ERS cuvette filled with co-substrate sam 0,1 mM
ERS cap filled with enzyme reagent with SAH Ado glutamate dehydrogenase

Stability and Storage

Kit can be used until expiry date printed on the labels, when stored at 2 - 8 °C. DO NOT FREEZE

Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Sample Material

Fresh serum or heparin plasma is recommended. Centrifuge samples immediately after collection to separate the serum from the blood cells. Haemolysed, turbid or severely lipaemic specimens are NOT recommended.

Reference Range

Adults: < 15 µmol/l
(Based on calibration against NIST SRM 1955 material)

It is recommended that every laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Homocysteine control kit should be used. Order number: ST1200

Calibration

Calibration has to be performed if quality control is out of range or the message *Calibration required* is shown on screen. For calibration the Eurolyser Homocysteine calibrator kit should be used. Order number: ST1210

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature before use! Test can be processed on instruments with **firmware version 1.36** or higher only.

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „measurement“ button, enter required information using the touch screen

2. Sample collection

- 2.1 Aspirate 20 µl sample from centrifuged tube.

3. Test processing

- 3.1 Dispense 20 µl sample INTO LIQUID of ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button.