hsCRP Testkit

Deutsch

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English

For human medical use only!

for quantitative in vitro determination of c-reactive protein in finger blood, venous blood, serum or plasma on the smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0105	hsCRP test kit	32 tests
Order number: ST1005	hsCRP control kit	1 x 1 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. A growing number of studies have examined that a more sensitive test, the so called hsCRP (high sensitive CRP) can predict recurrent cardiovascular disease, stroke and death in different settings. High levels of hsCRP consistently predict recurrent coronary events in patients with unstable angina and acute myocardial infarction (heart attack). Higher hsCRP levels also are associated with lower survival rates in these patients.

Method

Kinetic determination of the concentration of CRP by photometric measurement at 546 nm or 700 nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

Measurement Range

Immunoturbidimetric assay with 2 ranges:

Serum:	0,5 - 20 mg/l
Whole blood:	0,8 - 20 mg/l

Attention: LOT and HCT pending!

Test Kit

ERS cuvette pre-filled with buffer reagent

ERS cap pre-filled with latex reagent with rabbit anti human CRP antibodies

Stability and Storage

Storage at 2 - 8 °C

Unopened vacuum package: Until expiration date

3 months

Storage at 18 - 25 °C

Unrefrigerated: 1 month

Refrigerated at least 16h/day: 2 months

Warnings and Precautions

The reagents contain sodium azide (0,95 g/l) as preservative. DO NOT INGEST! Avoid contact with skin and eyes. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Whole blood or serum. Freshly drawn serum is preferred to guarantee highest test performance. *EDTA plasma samples show ~16% lower values, Li-Hep plasma samples show ~28% lower values compared to serum samples.*

Reference Range

Values < 1 mg/l	low cardiovascular risk
Values 1 - 3 mg/l	moderate cardiovascular risk
Values > 3 mg/l	elevated cardiovascular risk

It is recommended that each laboratory establishes its own reference range.

Quality Control

For internal quality control the hsCRP control kit should be used.

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ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

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- 1.** **Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Place ERS cap in test kit rack
 - 1.4 Press „Measurement“ button, enter required information using the touch screen

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- 2.** **Selection of whole blood or serum**

Please consult the laboratory photometer user manual for more information.

Sampletype menu:
Select the sample type
Haematocrit correction (if desired)
Automatic haematocrit-correction (**smart 700/546 and CUBE-S only!**)

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- 3.** **Processing for whole blood or serum**
 - 3.1a Aspirate 20 µl whole blood from fingertip using the minivette®
OR ...
 - 3.1b Aspirate 20 µl venous blood from EDTA primary tube
OR ...
 - 3.2 Aspirate 20 µl serum from primary tube
Dispense 20 µl sample INTO LIQUID in ERS cuvette
 - 3.3 Apply ERS cap firmly onto ERS cuvette
 - 3.4 Place ERS cartridge into laboratory photometer
 - 3.5 Start automatic sample processing by pressing the start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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