Lipoprotein (a) Control Kit

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Lipoprotein (a) Control Kit

for quality control procedures of the quantitative in vitro determination with the Lipoprotein (a) test kit on smart or CUBE (-S) laboratory photometer.

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

For human medical use only!

Order information

Order number: ST0141
Order number: ST1400

IVD

Indication

Lipoprotein (a) test kit
Lipoprotein (a) control kit

Kit size

16 Tests
2 x 1 ml (low/high)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This lipoprotein (a) control kit is intended for use in the clinical laboratory as a control for quantitative determination of lipoprotein (a) in the normal range and high range. Controls are prepared from a human serum and have been found to be non-reactive when tested with approved methods for HbsAg, anti HIV1/2 and anti HCV. However, the samples should be handled according to general rules in a clinical laboratory environment. This control kit is for in vitro diagnostic (IVD) use only.

Stability and Storage

The expiry date printed on the labels applies to storage of the unopened bottles at +2 to 8 °C

Stability after reconstitution:

Stored at +2 to 8 °C: 14 days

Preparation

Open the vial carefully and reconstitute the contents in 1 ml of distilled water by carefully rotating the vial (avoid frothing!). Allow reconstituted controls to stand for 10 minutes at room temperature before use. Treat the reconstituted control plasma like a citrated patient sample according to the instruction.

Waste Management

Please refer to your local legal requirements.

All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

Control values

Accuracy should be checked according to your local/national QC requirements with Lipoprotein (a) Control. The control plasma is treated in the same way as patient samples and results should be checked against the acceptable confidence limits supplied on the kit insert sheet and the vial label.

Warnings and Precautions

This testkit is for in vitro diagnostic use only. Do NOT INGEST. Avoid contact with skin and eyes. Take the necessary precautions for the use of laboratory reagents.


Measure the Controls same way like a patient sample and compare the values indicated in the lot-specific value sheet included in the control kit.

If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

Durchführung einer Lipoprotein (a) Kontrollmessung Processing of a Lipoprotein (a) control measurement

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
English

ATTENTION!

Allow single test to warm up at room temperature for a minimum of 10 minutes before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS R1 cuvette in the test kit rack
- 1.3 Place ERS R2 cap in the test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touch screen

2. Control preparation

- 2.1 Sex menu: Select „Control“

3. Processing control measurement with pipette

- 3.1 Aspirate 10 µl control liquid using a pipette
- 3.2 Dispense 10 µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS R2 cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer and start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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Value sheet for Lp(a) control kit

LOT

Niedrige Kontrolle - Low control

Ziel / Target	Wert / Value (min.)	Wert / Value (max)

Hohe Kontrolle - High control

Ziel / Target	Wert / Value (min.)	Wert / Value (max)