# Lipoprotein (a) test kit

**EURO**Lyser 

for quantitative in vitro determination of lipoprotein (a) on the smart or CUBE(-S) laboratory photometer.

 Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order information

Order number: ST0141  
Order number: ST1400

### Indication

Lipoprotein (a) test kit  
Lipoprotein (a) control kit

### Kit size

16 tests  
2 x 1 ml (low/high)

**Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

### Summary

The Eurolyser Lp(a) assay uses latex particles containing rabbit anti-human Lp(a) polyclonal antibody as a reagent. The rabbit anti-human Lp(a) polyclonal antibody is technically isoform sensitive by virtue of the antisera binding to multiple sites of kringle domain IV type 2 (KIV2); (The assay can theoretically be made nearly isoform independent if the appropriate calibrator system is used. This assay format, like most commercial assays, binds to both free apo(a) and true Lp(a) [i.e., apo(a) covalently bound to apoB-100]; therefore, it is best described as measuring "total apo(a)" rather than "Lp(a)."

The similarity of structural components of Lp(a) to LDL and to plasminogen suggests that Lp(a) may be associated with atherosclerosis and/or thrombosis. Although there is a lack of consistency in the conclusions of the studies about the contributory role of Lp(a) to coronary heart disease, it is widely accepted that Lp(a) is an important risk factor that may contribute to coronary artery disease independently or cooperatively with other risk factors. While the wide differences in Lp(a) levels seen among individuals are largely due to hereditary factors, the identification of individuals at risk through diagnostic screening should nevertheless be useful in alerting them to the need to eliminate or control other high risk factors when possible. Lp(a) values should be interpreted in conjunction with clinical evaluation and other lipoprotein tests when assessing atherosclerotic cardiovascular disease in specific populations.

### Method

Immunturbidimetric kinetic measurement of the Lp(a) concentration at 700 nm (Absorbance).

### Measurement Range

Measurement range: 0 - 100 mg/dl (LOT-pending)

Samples with concentrations higher than the upper limit of the measurement range must be diluted 1 + 2 with physiological saline (0.9% NaCl solution), e.g. 10 µl sample + 20 µl 0.9% NaCl solution, and the result multiplied by 3.

### Sample Material

Fresh sample material: Serum or EDTA plasma. Immediately centrifuge and process after blood collecting. Very lipaemic or turbid specimen must be clarified before assay is performed.

### Reference Range

Adults: < 30 mg/dl

It is recommended that every laboratory establishes its own reference ranges.

### Test Kit

ERS cuvette filled with glycine buffer solution  
ERS cap filled with latex reagent - rabbit anti human Lp(a) polyclonal antibodies

All reagent materials used are 100% original DENKA SEIKEN® raw materials.

### Quality Control

For internal quality control the Eurolyser Lipoprotein (a) control kit should be used. Order number: ST1400

### Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

### Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

### Waste Management

Please refer to local legal requirements.

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## English

### ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature before use!

#### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „measurement“ button, enter required information using the touch screen

#### 2. Sample collection

- 2.1 Aspirate 10 µl sample from centrifuged tube.

#### 3. Test processing

- 3.1 Dispense 10 µl sample INTO LIQUID of ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the smart laboratory photometer or by closing the door of the CUBE laboratory photometer.