# Microalbumin Testkit

## Deutsch

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# Microalbumin test kit

## English

for quantitative in vitro determination of albumin in urine on the smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0130	Microalbumin test kit	32 tests/pack
Order number: ST1300	Microalbumin control kit	2 x 1 ml (low/high)

 **Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

### Summary

A small amount of protein is excreted daily into the urine of healthy individuals. The extracted proteins are mucoproteins, most of which are filtered out of the uriniferous tubules and the glomeruli. Albumin, a protein of molecular weight of 50000 is not easily filtered out and is excreted into the urine (microalbuminuria) This makes albumin excretion into the urine a useful indicator of early glomerular disease. Microalbuminuria is a condition characterised by increased urinary excretion of albumin in the absence of overt nephropathy. Microalbumin has been reported in several studies to predict development of diabetic nephropathy and its mortality risk in diabetic patients. Because microalbuminuria may be reversible if diabetes is well controlled, the early detection of microalbumin may be very beneficial in treatment programs for diabetes.

### Principle

Antigen - Antibody reaction test; turbidity measurement at 700 nm.

### Measurement Range

Urine: 0.5 mg/dl - 30 mg/dl (5 mg/l - 300 mg/l)  
LOT dependent

### Test Kit

ERS cuvette filled with buffer reagent pH 7.6 100 mM Tris  
ERS cap filled with antiserum reagent pH 7.6  
20% anti-human albumin goat antiserum 100 mM Tris

### Storage and Stability

Store at 2 - 8 °C. Kit can be used until expiry date printed on the labels.  
DO NOT FREEZE!

### Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST! Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

### Waste Management

Please refer to local legal requirements.

### Sample Material

The specimen should be fresh or 24 hour urine. Urine specimen should be stored refrigerated 2 - 8 °C. The specimen may be stored refrigerated up to 2 weeks.

### Reference Range

The expected value for microalbumin is < 2.5 mg/dl (< 25 mg/l). It is recommended that each laboratory establishes its own reference ranges.

### Quality Control

For internal quality control the Microalbumin control kit should be used.  
Order number: ST1300

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### ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature before use!

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#### 2. Test preparation

- 2.1 Aspirate 20 µl sample (urine) out of sample tube with pipette
- 2.2 Dispense 20 µl sample into ERS cuvette INTO LIQUID

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#### 3. Test processing

- 3.1 Apply ERS cap firmly onto ERS cuvette
- 3.2 Place ERS cartridge into laboratory photometer
- 3.3 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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