bovine Progesterone

VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of bovine Progesterone in li-hep plasma or serum on a solo or CUBE-VET analyser

 Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order information

Order number: VT0260
Order number: VT2600

Indication

bovine Progesterone VET test kit
bovine Progesterone VET control kit

Kit size

6 tests/pack
1 x 2 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The progesterone level changes during the sexual cycle of cows. At day 0 of the cycle no or only very low levels of progesterone can be detected due to an inactive corpus luteum. The progesterone levels reach their maximum around 13 – 14 days after the estrus. In case the cow is pregnant, progesterone levels increase continuously up to day 21 after fertilisation. Decrease of progesterone levels at day 20 show the begin of a new cycle which can be used for a new insemination. Therefore, measuring of progesterone levels is a valuable tool for cattle farms to improve the reproductive management.

Method

Homogeneous immunoturbidimetric test.

Measurement Range

1.4 – 10.0 ng/ml

Sample Material

Use 40 µl of li-hep plasma or serum.

- Use fresh sample material and perform the measurement within 2 hours.
- Storage up to 24 hours: Separate serum/plasma from the cells and store it in a plain tube at 2 – 8 °C.
- Storage up to 2 weeks: Separate serum/plasma from the cells and store it in a plain tube at -20 °C

Test Kit

ERS cuvette filled with antibody reagent
ERS cap filled with buffer reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Decision Limits

Cow: > 2 ng/ml (gestation)

It is recommended that each laboratory establishes its own decision limits.

Quality Control

For internal quality control the Eurolyser bovine Progesterone VET control kit is recommended. Order number: VT2600

Precision

Reproducibility within-run:

Control; N = 20; mean = 5.5 ng/ml; CV = 8.51%

Precision accepted: CV < 12%

Correlation

Bovine sample correlation: N = 50

y (Eurolyser bProg) = 1.0005x (Siemens Progesterone) - 0.0124;
 $R^2 = 0.9424$;

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 15% of initial values

Haemoglobin	525 mg/dl
Human albumin	12 g/dl
Bilirubin (conjugated)	72 mg/dl
Bilirubin (unconjugated)	72 mg/dl
Cholesterol	620 mg/dl
Rheumatoid factor	1080 IU/ml
Triglycerides	835 mg/dl
Uric acid	30 mg/dl

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touchscreen

2. Sample preparation

- 2.1 Aspirate 40 µl sample material from centrifuged sample tube

3. Sample processing

- 3.1 Dispense 40 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.