Cortisol

VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of cortisol in li-hep plasma or serum on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
 Bindergasse 3
 5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
 Fax +43 (0)662 / 43 21 00-50
 www.eurolyser.com

English

Order information

Order number: VT0290
 Order number: VT0291
 Order number: VT2900

Indication

Cortisol VET test kit
 Cortisol VET test kit
 Cortisol VET control kit

Kit size

16 tests
 6 tests
 1 x 0.5 ml (decision level)

Important: AFTER preparing the reagent, allow single test at least 10 min. to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Cortisol is a glucocorticoid hormone that affects the metabolism of carbohydrates, proteins, and fats. Chronic, excessive glucocorticoid release leads to hyperadrenocorticism (HAC). Pituitary-dependent hyperadrenocorticism (PDH), more commonly known as Cushing's disease, is caused by a pituitary tumour (mostly adenomas) in the brain. The tumour triggers excessive levels of adrenocorticotrophic hormone (ACTH), which stimulates the adrenal glands near the kidneys to produce cortisol. A small percentage of dogs with HAC have a tumour of one of the adrenal glands. This form of HAC is called adrenal-dependent HAC and results from a direct increase in cortisol production by the adrenal gland tumour. The overproduction of cortisol causes symptoms such as hair loss, pot-bellied appearance, increased appetite, and polydipsia and polyuria. Addison's disease (hypoadrenocorticism, or underproduction of cortisol) occurs less commonly than the Cushing's disease in dogs.

Method

Homogeneous, immunoturbidimetric test

Measurement Range

0.8 – 24 µg/dl

Sample Material

Use 40 µl of li-hep plasma or serum

- Use fresh sample material and perform the measurement within 2 hours.
- Storage up to 24 hours: Separate serum/plasma from the cells and store it in a plain tube at 2 – 8 °C.
- Storage up to 2 weeks: Separate serum/plasma from the cells and store it in a plain tube at -20 °C

Test Kit

ERS cuvette filled with antibody reagent
 ERS cap filled with buffer reagent
 Bottle filled with dilution buffer

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 15% of initial values.

Haemoglobin	131	mg/dl
Albumin	12	g/dl
Bilirubin, conjugated	36	mg/dl
Bilirubin, unconjugated	36	mg/dl
Cholesterol	620	mg/dl
Rheumatoid factor	540	IU/ml
Triglycerides	835	mg/dl
Uric acid	30	mg/dl

Precision

Reproducibility within-run:
 Control; N = 20; mean = 5.7 µg/dl; CV = 6.93%

Precision accepted: CV < 15%

Correlation

Sample correlation: N = 47

y (Eurolyser cortisol) = 0.9168x (Roche Elecsys) - 0.2049; R² = 0.9624;

Quality Control

For internal quality control the Eurolyser Cortisol VET control kit is recommended. Order number: VT2900

Reference Ranges

It is recommended that each laboratory establishes its own reference ranges.

I. ACTH Stimulation Test

Pre-ACTH	Post-ACTH	Interpretation
< 2 µg/dl	< 2 µg/dl	Consistent with hypoadrenocorticism
	2 – 6 µg/dl	Inconclusive
2 – 6 µg/dl	6 – 18 µg/dl	Normal
	18 – 22 µg/dl	Equivocal, Cushing's syndrome possible
> 22 µg/dl	> 22 µg/dl	Consistent with Cushing's syndrome

II. Low-dose Dexamethasone Suppression Test

4-hour cortisol level	8-hour cortisol level	Interpretation
< 1 µg/dl	< 1 µg/dl	Normal
1 – 1.5 µg/dl	1 – 1.5 µg/dl	Inconclusive
> 1.5 µg/dl and > 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with Cushing's syndrome
< 1.5 µg/dl or < 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with PDH
> 1.5 µg/dl or > 50% of baseline	> 1.5 µg/dl and < 50% of baseline	Consistent with PDH

III. High-dose Dexamethasone Suppression Test

4-hour cortisol level	8-hour cortisol level	Interpretation
< 1.5 µg/dl or < 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with PDH
> 1.5 µg/dl and > 50% of baseline	< 1.5 µg/dl or < 50% of baseline	Consistent with PDH
< 1.5 µg/dl or < 50% of baseline	< 1.5 µg/dl or < 50% of baseline	Consistent with PDH
> 1.5 µg/dl and > 50% of baseline	> 1.5 µg/dl and > 50% of baseline	Consistent with PDH or ATH

IV. Post ACTH Stimulation Mitotane (Lysodren®) Dosing and Monitoring

Loading Dose	Interpretation	Maintenance Dose	Interpretation
< 1 µg/dl	Discontinue mitotane	< 1 µg/dl	Discontinue mitotane
1 – 5 µg/dl	Begin maintenance mitotane dosing	1 – 5 µg/dl	Maintain current dosage
> 5 µg/dl	Continue mitotane loading dose	> 5 µg/dl	Increase weekly dose or repeat loading dose.

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English

ATTENTION!

Tests can be processed on **solo** laboratory photometers with version **2.15** (or higher) or **CUBE-VET** laboratory photometers with version **4.1** (or higher)!

AFTER preparing the reagent, allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Press „Measurement“ button, enter required information using the touchscreen
- 1.3 Place ERS cuvette in test kit rack
- 1.4 Place ERS cap in test kit rack
- 1.5 Take out bottle containing dilution buffer from the test kit package

Allow single test to warm up to room temperature for at least 10 minutes before measurement! Put test kit package back into refrigerator

- 2.3 Aspirate 40 µl sample material from centrifuged sample tube

3. Sample processing

- 3.1 Dispense 40 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.