canine Progesterone

VET test kit

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of canine Progesterone in li-hep plasma or serum on a solo or CUBE-VET analyser

English

	Order information	Indication	Kit size
	Order number: VT0230 Order number: VT0231 Order number: VT2300	canine Progesterone VET test kit canine Progesterone VET test kit canine Progesterone VET control kit	6 tests 16 tests 1 x 2 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Progesterone (P4) is an androgen produced by the corpus luteum of bitches. Normally present in low doses in non-pregnant bitches, in dogs, levels of progesterone start to climb around the time of the LH peak. Thus, the production of progesterone by the follicle begins prior to ovulation at the end of proestrus and reaches increased levels in peripheral blood plasma at the time of ovulation. In females, the measurement of Progesterone is useful in evaluating the status of ovarian functions, monitoring of Progesterone therapy and early stage pregnancy evaluations.

Method

Homogeneous immunoturbidimetric test.

Measurement Range

When using a 40 µl pipette: 1.2 – 8.0 ng/ml
(3.8 – 25.4 nmol/l)

When using a 5 µl pipette: 6.0 – 20.0 ng/ml
(19.1 – 63.6 nmol/l)

Sample Material

Li-hep plasma or serum. Sample volume according to the menu settings on the laboratory photometer and according to the use of a 40 µl or 5 µl pipette.

- Use fresh sample material and perform the measurement within 2 hours.
- Storage up to 24 hours: Separate serum/plasma from the cells and store it in a plain tube at 2 – 8 °C.
- Storage up to 2 weeks: Separate serum/plasma from the cells and store it in a plain tube at -20 °C

Test Kit

ERS cuvette filled with antibody reagent
ERS cap filled with buffer reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Decision Limits

Dog: 4 – 8 ng/ml (12.7 – 25.4 nmol/l) (Ovulation)

It is recommended that each laboratory establishes its own decision limits.

Quality Control

For internal quality control the Eurolyser canine Progesterone VET control kit is recommended. Order number: VT2300

Precision

Reproducibility within-run:
Control; N = 20; mean = 3.5 ng/ml; CV = 9.57%

Precision accepted: CV < 12%

Correlation

Canine sample correlation: N = 50

y (Eurolyser progesterone) = 1.0098x (Immulite progesterone) - 0.0028; $R^2 = 0.945$;

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 15% of initial values

Haemoglobin	525 mg/dl
Human albumin	12 g/dl
Bilirubin (conjugated)	72 mg/dl
Bilirubin (unconjugated)	72 mg/dl
Cholesterol	620 mg/dl
Rheumatoid factor	1080 IU/ml
Triglycerides	835 mg/dl
Uric acid	30 mg/dl

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touchscreen

2. Sample preparation

Choose sample type:

- 2.1a „Pipette“ for measurement range:
1.2 - 8.0 ng/ml
OR ...
- 2.1b „Pipette5µ“ for measurement range:
6.0 - 20.0 ng/ml
- 2.2 Aspirate 40 µl **OR** 5 µl sample material from centrifuged sample tube

3. Sample processing

- 3.1 Dispense sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.