equine Progesterone

VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of equine Progesterone in serum or lithium heparin plasma on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order information

Order number: VT 0250
Order number: VT 2500

Indication

equine Progesterone VET test kit
equine Progesterone VET control kit

Kit size

6 tests/pack
1 x 2 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

During early pregnancy, progesterone is produced in the equine ovary by the corpus luteum (CL). Its concentration remains elevated and peaks between 60 and 120 days of gestation. Circulating progesterone is used diagnostically to evaluate luteal function during early pregnancy. When the circulation progesterone concentration reaches a distinct level, this can be considered adequate to maintain early stage pregnancy. There are a number of reasons for monitoring during pregnancy, such as uterine infections, history of pregnancy loss and luteal insufficiency.

Method

Homogeneous immunoturbidimetric test.

Measurement Range

1.4 – 13.0 ng/ml

Sample Material

Use 40 µl of serum or li-hep plasma.

Test Kit

ERS cuvette filled with buffer reagent
ERS cap filled with antibody reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Decision Limits

Horse: > 2 ng/ml (gestation)

It is recommended that each laboratory establishes its own decision limits.

Quality Control

For internal quality control the Eurolyser equine Progesterone VET control kit is recommended. Order number: VT 2500

Precision

Reproducibility within-run:

Control; N = 20; mean = 5.4 ng/ml; CV = 9.69%

Correlation

Equine sample correlation: N = 50

y (Eurolyser eProg) = 1.0055x (Siemens Progesterone) - 0.0119;
 $R^2 = 0.9522$;

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 15% of initial values

Haemoglobin	525 mg/dl
Human albumin	12 g/dl
Bilirubin (conjugated)	72 mg/dl
Bilirubin (unconjugated)	72 mg/dl
Cholesterol	620 mg/dl
Rheumatoid factor	1080 IU/ml
Triglycerides	835 mg/dl
Uric acid	30 mg/dl

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


Processing of a equine Progesterone VET test

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





English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touchscreen

2. Sample preparation

- 2.1 Aspirate 40 µl sample material from centrifuged sample tube

3. Sample processing

- 3.1 Dispense 40 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.