# HGB VET control kit

**For veterinary use only!**

For quality control procedures of the quantitative determination with the HGB VET test kit on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order information

Order number: VT 0150  
Order number: VT 1500

### Indication

HGB VET test kit  
HGB VET control kit

### Kit size

32 tests/pack  
1 x 1.5 ml (decision level)

**Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

This HGB control kit is intended for use in the clinical laboratory as a control for quantitative determination of HGB. It is prepared from haemoglobin equivalent.

## Preparation and Handling

The controls are liquid and ready-to-use. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. We recommend not to pipette directly from the original bottle as well as using single use pipetting tips to avoid contamination. Dispense about 5 drops into a serum cup and pipette out of it. Dispose cup after use. For further details on how to use the control kit please refer to the description inside this package insert. Turbidity shows bacterial contamination, the vial has to be discarded.

## Stability and Storage

Unopened vials are stable at 4 - 25 °C until the expiration date stated on the label. Once opened vials are stable until the expiry date when stored tightly capped at 4 - 25 °C.

## Waste Management

Please refer to local legal requirements

## Warnings and Precautions

This control kit is for veterinary diagnostic use only!  
**DO NOT INGEST!** Exercise the normal precautions required for handling all laboratory reagents. All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient specimen. In the event of exposure the directives of the responsible health authorities should be followed.

Set the species to „Control“ before running the QC material

Compare the results with the values indicated in the lot-specific value table

If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor

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


## Processing of a HGB VET control measurement

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## English

### ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

#### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Detach cap from ERS cuvette
- 1.4 Press „Measurement“ button, enter required information using the touchscreen

#### 2. Control preparation

*For further details please see user manual of laboratory photometer*

- 2.1 Set species to „Control“
- 2.2 Aspirate 20 µl control liquid using a pipette

#### 3. Control processing

- 3.1 Dispense 20 µl control liquid INTO THE LIQUID in the ERS cuvette
- 3.2 Apply cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

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Technical details subject to change without notice.

Doc-ID: fcdfb139-9853-4ab4-9bf7-546ce4fe2026 ; latest revision: DP 1.0 ; date of print: 2018-09-13