HGB VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of haemoglobin in whole blood on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order information

Order number: VT 0150
Order number: VT 1500

Indication

HGB VET test kit
HGB VET control kit

Kit size

32 tests/pack
1 x 1 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

In routine medical diagnostics, the test may be used for the following purposes: To diagnose anaemia or polyglobulia, to monitor risk groups for iron deficiency, as part of point-of-care testing when acute bleedings or haemorrhages are suspected, at parasitic diseases as well as for quantifying blood loss.

Method

The reagent binds haemoglobin and its variants to form a coloured complex, which is detected photometrically at 546 nm.

Measurement Range

Measurement range: 0 - 30 g/dl
Lower detection limit: 0.5 g/dl

Sample Material

20 µl whole blood in anticoagulant tube (EDTA, heparin).
Make sure sample is mixed well.

Test Kit

Filled with AHD reagent.

Stability and Storage

Store at 2 - 25 °C. Kit can be used until expiry date printed on the labels. DO NOT FREEZE!

Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST! Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Reference Range

Dog	15.0 - 20.0 g/dl	(9.30 - 12.40 mmol/l)
Cat	8.0 - 17.0 g/dl	(4.97 - 10.56 mmol/l)
Horse	9.0 - 18.0 g/dl	(5.59 - 11.18 mmol/l)
Rabbit	11.5 - 15.1 g/dl	(7.14 - 9.38 mmol/l)
Pig	11.0 - 16.5 g/dl	(6.83 - 10.25 mmol/l)
Sheep	11.0 - 16.5 g/dl	(6.83 - 10.25 mmol/l)
Cattle	10.0 - 13.0 g/dl	(6.21 - 8.07 mmol/l)
Mouse	11.1 - 18.0 g/dl	(6.89 - 11.18 mmol/l)
Rat	11.1 - 18.0 g/dl	(6.89 - 11.18 mmol/l)

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser HGB VET control kit is recommended. Order number: VT 1500

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Detach cap from ERS cuvette
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

-
- 2. Sample preparation**
 - 2.1 Aspirate 20 µl sample material from EDTA primary tube

-
- 3. Sample processing**
 - 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.