K⁺ Potassium VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of K⁺ Potassium in serum or lithium-heparin plasma on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT 0280 Order number: VT 0281 Order number: VT 2800	K ⁺ Potassium VET test kit K ⁺ Potassium VET test kit K ⁺ Potassium VET control kit	16 tests/kit 6 tests/kit 1 x 5 ml (decision level)

Summary

Potassium is the most important cation in the intracellular liquid. It represents a crucial buffer within the cell. Potassium facilitates nerve conduction as well as muscle function and maintains the osmotic pressure. Abnormal potassium levels influence the breathing muscle stimulation and heart muscle function.

Elevated potassium levels (hyperkalaemia) occur with kidney insufficiency, urinary tract obstruction, respiratory or metabolic acidosis and hypoadrenocorticism, as well as excessive haemolysis in cats, horses, cattle and several dog breeds.

Usually, reduced potassium levels (hypokalaemia) are the result of excessive salt loss caused by severe diarrhea or regurgitation.

Low potassium values on the other hand correlate with insufficient food intake (anorexia, especially in cats), malabsorption and severe burns.

Method

Immunturbidimetric test: The reagent agglutinates potassium ions and forms a stable, cloudy suspension. The degree of turbidity is measured by the photometer at 546 nm.

Measurement Range

0.6 - 12 mmol/l

Sample Material

20 µl of fresh serum (not older than 1 hour) or lithium-heparin plasma can be used. Results from li-hep plasma samples are about 8% lower than results from serum samples. Haemolytic or lipaemic samples must never be used!

Test Kit

ERS cuvette prefilled with sodium hydroxide / sodium teraphenylborate and stabilisers.

Stability and Storage

Stable until the expiration date stated on the label when stored at 15 - 25 °C.

DO NOT STORE COOLED AND DO NOT FREEZE!

Warnings and Precautions

This test kit is for in vitro diagnostic use only! Do not swallow! Avoid contact with skin and eyes. Reagents contains sodium azide as preservative (0.95 g/l), which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Reference Range

Dog:	3.5 - 5.1 mmol/l
Cat:	3.0 - 4.8 mmol/l
Rabbit:	3.7 - 6.3 mmol/l
Guinea pig:	4.5 - 8.8 mmol/l
Cattle:	3.5 - 4.5 mmol/l
Horse:	2.8 - 4.5 mmol/l
Pig:	4.0 - 5.0 mmol/l

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser K⁺ Potassium VET control kit is recommended. Order number: VT 2800

Precision

Reproducibility within-run:

Control 1: N = 20; mean = 4.87 mmol/l; CV = 3.11%

Control 2: N = 20; mean = 7.75 mmol/l; CV = 4.20%

Correlation

N = 38; $y(\text{Eurolyser K}^+) = 1.1772x(\text{Spotchem EL}) - 1.0489$;
 $R^2 = 0.9459$

Interferences

Haemolytic and lipaemic samples produce incorrectly elevated K⁺ values.

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Bezeichnung

K+ Potassium VET Testkit
K+ Potassium VET Testkit
K+ Potassium VET Control Kit

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16 Tests/Packung
6 Tests/Packung
1 x 5 ml (decision level)

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


Processing of an K⁺ Potassium VET test

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








English

ATTENTION!

Tests can be processed on **solo** laboratory photometers with the following version:

- Black and white display: **1.42** (or higher)
- Colour display: **2.09** (or higher)

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Instrument ready to start the reagent blank measurement.

2. Reagent Blank Measurement

- 2.1 Place ERS cuvette WITH cuvette cap and WITHOUT sample into the laboratory photometer. Start automatic reagent blank measurement by pressing the  start button on the smart or by closing the door of the CUBE-VET laboratory photometer.
- 2.2 After the reagent blank measurement take the ERS cartridge out of the laboratory photometer and confirm by pressing the  /  button.
- 2.3 Instrument ready for start of analysis.

3. Sample processing

- 3.1 Detach cap from ERS cuvette
- 3.2 Aspirate 20 µl sample material from centrifuged sample tube
- 3.3 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
- 3.4 Apply cap firmly onto ERS cuvette
- 3.5 Place ERS cartridge into laboratory photometer. Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.