Lactate VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of lactate in li-hep plasma on a solo or CUBE-VET analyser.

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT0220 Order number: VT0221 Order number: VT2200	Lactate VET test kit Lactate VET test kit Lactate VET control kit	16 tests/pack 6 tests/pack 1 x 5 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

In the face of anaerobic metabolism lactate is synthesized from pyruvate. Muscle tissue is the predominant source of lactate production. Increased lactate concentration is a sensitive indication for a significant metabolic disturbance. Possible reasons are hypoxia and hypoperfusion. Diseases associated with increased lactate concentration are hypovolemic shock, local hypoperfusion (i.e. intestinal ischaemia in the face of ileus, aortic-thromboembolism), hypoxia, anaemia, or sepsis. In racehorses lactate is used to either check performance ability or i.e. exclude latent myopathy.

Method

Colorimetric test. Photometric endpoint measurement of the absorbance at 546 nm wavelength.

Measurement Range

0.1 – 15 mmol/l

When values exceed this range, samples should be diluted 1 + 1 with NaCl solution (9 g/l) and the results multiplied by 2.

Sample Material

Use 20 µl of li-hep plasma. Do not use serum!

For measurements right after taking the blood sample:

Use lithium heparin samples and separate plasma from red blood cells immediately (within 5 minutes).

In case measurements cannot be done immediately:

Use sodium fluoride/EDTA samples, centrifuge and measure at the earliest.

Stability in plasma (after centrifugation and separation):

At 15 – 25 °C: 8 hours

At 2 – 8 °C: 14 days

Test Kit

Filled with mono reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Reference Range

Canine:	0.5 – 3.0 mmol/l
Feline:	< 1.0 mmol/l
Equine:	0.5 – 2.0 mmol/l

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Lactate VET control kit is recommended. Order number: VT2200

Precision

Reproducibility within-run:

Control high; N = 20; mean = 3.87 mmol/l; CV = 1.69%;

Control low; N = 20; mean = 1.49 mmol/l; CV = 2.10%;

Precision accepted: CV < 10%

Correlation

Equine sample correlation: N = 40

y (Eurolyser Lactate) = 1.0602 x (cobas Lactate) - 0.2564; R^2 = 0.9577;

Interferences

No interferences up to:

Ascorbic acid	30 mg/dl
Bilirubin	60 mg/dl
Triglyceride	2000 mg/dl
Haemoglobin	1000 mg/dl
Dopamine	10 mg/l
L-Dopamine	20 mg/l
Methyldopamine	10 mg/l
Glycolic acid	1200 mg/l

For further information on interfering substances refer to ¹

References

1. YOUNG, D.S. (2000) Effects of drugs on clinical laboratory tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press;

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Detach cap from ERS cuvette
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

-
- 2. Sample preparation**
 - 2.1 Aspirate 20 µl sample material from centrifuged sample tube

-
- 3. Sample processing**
 - 3.1 Dispense 20 µl sample INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply ERS cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.