NH₃ Ammonia VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of NH₃ Ammonia in li-hep plasma or EDTA plasma on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order information

Order number: VT0160
Order number: VT1600

Indication

NH₃ Ammonia VET test kit
NH₃ Ammonia VET control kit

Kit size

6 tests/pack
1 x 5 ml

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Protein digestion produces ammonia within the intestine by deamination of amino acid. Ammonia reaches the liver by portal blood where it is synthesized into urea (main route of metabolism of intestinal ammonia into urea).

Heavy liver dysfunction (e.g. shunt, cirrhosis) can prevent the liver from synthesizing ammonia into urea. As a result one will receive a lowered urea concentration in relation to creatinine as well as increased blood ammonia. Ammonia and other protein metabolism (aromatic amino acids) that get out of control of the liver metabolites can affect the central nervous system (hepato-encephalic syndrome).

Method

Enzymatic test. Photometric measurement of the absorbance rate (kinetic) at 340 nm wavelength.

Measurement Range

10 – 1100 µmol/l (Lot dependent)

Sample Material

Use only 40 µl fresh li-hep plasma or EDTA plasma (Ammonium heparin plasma or serum must NOT be used)! After drawing the blood, immediately centrifuge and analyse within 5 minutes on the analyser.

Test Kit

Filled with buffer, pH 8.0, alpha-ketoglutarate, GLDH and NADPH.

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C.

Opening the vacuum package may limit the reagent stability to three months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Reference Range

Dog/Cat: < 60 µmol/l

It is recommended that each laboratory establishes its own reference range.

Quality Control

For internal quality control the Eurolyser NH₃ Ammonia VET control kit is recommended. Order number: VT1600

Precision

Reproducibility within-run:

Control; N = 20; mean = 245 µmol/l; CV = 4.89%;

Precision accepted: CV < 12%

Correlation

y (Eurolyser NH₃) = 0.9731 x (Pentra NH₃) + 0.5393;
R² = 0.9759

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Ascorbic acid	3 mg/dl
Bilirubin	40 mg/dl
Lactate	200 mg/dl
Pyruvate	6.6 mg/l
Triglycerides	200 mg/dl

References

1. HARDY, RM. (1989) Diseases of the liver and their treatment, In: Ettinger SJ, ed. Textbook of veterinary internal medicine, 3rd ed. Philadelphia: WB Saunders; 1479 – 1527;
2. RASKIN, NH., PRICE, JB., FISHMAN, RA. (1964) Portal-systemic encephalopathy due to congenital intrahepatic shunts. N Engl J Med; 270:225 – 229;
3. RASKIN, NH., BREDESEN, D., EHRENFIELD, WK., et al. (1984) Periodic confusion caused by congenital extrahepatic portacaval shunt; Neurology; 34:666 – 669;

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

Referenzen

1. HARDY, RM. (1989) Diseases of the liver and their treatment, In: EttingerSJ, ed. Textbook of veterinary internal medicine, 3rd ed. Philadelphia: WB Saunders; 1479 – 1527;
2. RASKIN, NH., PRICE, JB., FISHMAN, RA. (1964) Portal-systemic encephalopathy due to congenital intrahepatic shunts. N Engl J Med; 270:225 – 229;
3. RASKIN, NH., BREDESEN, D., EHRENFIELD, WK., et al. (1984) Periodic confusion caused by congenital extrahepatic portacaval shunt;

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
Processing of an NH₃ Ammonia VET test

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

- 1. Preparation of test system**
 - 1.1 Place RFID card
 - 1.2 Place ERS cuvette in test kit rack
 - 1.3 Detach cap from ERS cuvette
 - 1.4 Press „Measurement“ button, enter required information using the touchscreen

- 2. Sample preparation**
 - 2.1 Aspirate 40 µl sample material from centrifuged sample tube

- 3. Sample processing**
 - 3.1 Dispense 40 µl sample INTO THE LIQUID in the ERS cuvette
 - 3.2 Apply ERS cap firmly onto ERS cuvette
 - 3.3 Place ERS cartridge into laboratory photometer
 - 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.