Total Bile Acids

VET test kit For veterinary use only!

EUROLyser

Intended use:

Veterinary test kit for quantitative in vitro determination of bile acids in serum or lithium heparin plasma on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

Order information

Order number: VT0240
Order number: VT0241
Order number: VT2400

Indication

Total Bile Acids VET test kit
Total Bile Acids VET test kit
Total Bile Acids VET control kit

Kit size

16 tests
6 tests
1 x 2 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Bile acids are an important parameter to evaluate liver function. Increased bile acid concentrations can occur due to various hepatopathies, portosystemic shunts or intra- or extrahepatic cholestasis. The bile acid concentration depends on several factors (i.e. food ingestion, gall bladder emptying). Therefore, it is recommended to only analyse bile acids if the patient is fasting or after performance of a bile acid stimulation test.

Method

Homogeneous, colorimetric test.

Measurement Range

3 – 200 µmol/l

Sample Material

Use 100 µl of serum or li-hep plasma.
Lypaemic and haemolytic samples shall not be used.

Test Kit

ERS cuvette filled with colour reagent
ERS cap filled with activator reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to 3 months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes.
Observe all necessary precautions for the use of laboratory reagents.

Single determinations can lead to false positive results, therefore the performance of bile acid stimulation tests is preferable, except in horses.

Please note that animals must fast 12 hours before sampling (except horses).

Waste Management

Please refer to local legal requirements.

Reference Ranges

Dog: < 20 µmol/l
Cat: < 20 µmol/l
Horse: < 12 µmol/l

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser Total Bile Acids VET control kit is recommended. Order number: VT 2400

Precision

Control; N = 20; mean = 29.12 µmol/l; CV = 2.57%

Correlation

Sample correlation: N = 72

y (Eurolyser TBA) = 0.9924x (reference TBA) + 0.2924;
 $R^2 = 0.9774$;

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Bilirubin:	35 mg/dl
Haemoglobin:	120 mg/dl
Triglycerides:	1200 mg/dl

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

Processing of a Total Bile Acids VET test

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





English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touchscreen

2. Sample preparation

- 2.1 Aspirate 100 µl sample material from centrifuged sample tube

3. Sample processing

- 3.1 Dispense 100 µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.