# GLDH VET test kit

## For veterinary use only!

Veterinary test kit for quantitative in vitro determination of GLDH in li-hep plasma or serum on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH  
Bindergasse 3  
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00  
Fax +43 (0)662 / 43 21 00-50  
www.eurolyser.com

## English

### Order information

Order number: VT0200  
Order number: VT0201  
Order number: VT2000

### Indication

GLDH VET test kit  
GLDH VET test kit  
GLDH VET control kit

### Kit size

16 tests/pack  
6 tests/pack  
1 x 5 ml (decision level)

**Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

Glutamate dehydrogenase (GLDH) is a liver-specific enzyme which is located in mitochondria of centrolobular hepatocytes. Due to this localization, GLDH dominantly increases in hypoxic states of the liver as well as in pronounced necrosis of hepatocytes.

Increased release of GLDH from hepatocytes can further be observed in i.e. inflammatory, neoplastic, or infectious liver pathology as well as in all diseases leading to liver cell degeneration.

## Method

GLDH catalyses the reduction of oxo-glutarate and at the same time the oxidation of NADH to NAD.

The reaction speed of the absorption increase is directly proportional to the GLDH activity which is measured at a wavelength of 340 nm.

## Measurement Range

6 – 120 U/L (Lot dependent)

## Sample Material

Use 100 µl li-hep plasma or serum samples.  
(Avoid lipaemic samples)

In case samples need to be kept, store in a dark place only!

### Sample material stability:

Stored at 20 °C:	4 hours
Stored at 2 – 8 °C:	1 day
Stored at -20 °C:	4 weeks

## Test Kit

ERS cuvette filled with: Triethanolamine, alpha-ketoglutarate, ammonium acetate, EDTA, ADP, LDH.  
ERS cap filled with: NADH

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

## Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

## Waste Management

Please refer to local legal requirements.

## Reference Range

Dog:	< 13 U/L
Cat:	< 10 U/L
Rabbit:	< 19 U/L
Guinea pig:	< 17 U/L
Horse:	< 13 U/L
Bovine:	< 15 U/L
Sheep:	< 11 U/L
Pig:	< 7 U/L
Goat:	< 20 U/L

It is recommended that each laboratory establishes its own reference ranges.

## Quality Control

For internal quality control the Eurolyser GLDH VET control kit is recommended. Order number: VT2000

## Precision

Reproducibility day-to-day:  
Control; N = 20; mean = 24.4 U/L; CV = 8.02%;

Precision accepted: CV < 15%

## Correlation

y (Eurolyser GLDH) = 1.0414 x (Cobas GLDH) + 4.2329; R<sup>2</sup> = 0.9494.

## Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Ascorbic acid	30 mg/dl
Bilirubin	60 mg/dl
Haemoglobin	500 mg/dl
Triglyceride	400 mg/dl

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### Bezeichnung

GLDH VET Testkit  
GLDH VET Testkit  
GLDH VET Control Kit

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16 Tests/Packung  
6 Tests/Packung  
1 x 5 ml (decision level)

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---

### ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) before use!

- 1. Preparation of test system**
  - 1.1 Place RFID card
  - 1.2 Place ERS cuvette in test kit rack
  - 1.3 Place ERS cap in test kit rack
  - 1.4 Press „Measurement“ button, enter required information using the touchscreen

- 
- 2. Sample preparation**
    - 2.1 Aspirate 100 µl sample material from centrifuged sample tube

- 
- 3. Sample processing**
    - 3.1 Dispense 100 µl sample INTO THE LIQUID in the ERS cuvette
    - 3.2 Apply ERS cap firmly onto ERS cuvette
    - 3.3 Place ERS cartridge into laboratory photometer
    - 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.