

Human test information sheet

PARAMETER	SAMPLE MATERIAL	SAMPLE VOLUME
CRP	Finger blood / EDTA venous blood / Serum	5 µl
hsCRP	Finger blood / EDTA venous blood / Serum	20 µl
Cystatin C (GFR)	Finger blood / EDTA venous blood / Serum	20 µl
D-dimer	Citrat / EDTA plasma	20 µl
iFOB (faecal occult blood)	Faeces	40 µl
Ferritin	Serum / Li-hep plasma from capillary tube	50 µl
HbA1c	Finger blood / EDTA venous blood	10 µl
Lipoprotein (a)	Serum / EDTA plasma	10 µl
PT (INR)	Finger blood / Citrat venous blood / Citrat plasma	20 µl
Troponin I	Serum / Li-hep plasma	50 µl