

Human test information sheet

PARAMETER	SAMPLE MATERIAL	SAMPLE VOLUME	CUBE-S	CUBE	smart 700/546	smart 700/340
ASO	Finger blood / EDTA venous blood / Serum	5 µl	●	●	●	●
CRP	Finger blood / EDTA venous blood / Serum	5 µl	●	●	●	●
hsCRP	Finger blood / EDTA venous blood / Serum	20 µl	●	●	●	●
Cystatin C	Finger blood / EDTA venous blood / Serum	20 µl	●	●	●	●
D-dimer	Citrat / EDTA plasma	20 µl	●	●	●	●
iFOB (faecal occult blood)	Faeces	40 µl	●	●	●	●
Ferritin	Serum / Li-hep plasma from capillary tube	50 µl	●	●	●	●
Haemoglobin	Finger blood / EDTA venous blood	20 µl	●		●	
Homocysteine	Serum / Li-hep plasma	20 µl				●
HbA1c	Finger blood / EDTA venous blood	10 µl	●	●	●	●
K+ Potassium	Serum / Li-hep plasma	20 µl	●	●	●	●
Lipoprotein (a)	Serum / EDTA plasma	10 µl	●	●	●	●
NH3 Ammonia	EDTA plasma	40 µl				●
Microalbumin	Urine	20 µl	●	●	●	●
PT (INR)	Finger blood / Citrat venous blood / Citrat plasma	20 µl	●	●	●	●
Troponin I	Serum / Li-hep plasma	50 µl	●	●	●	●