

VET test information sheet

PARAMETER	SAMPLE MATERIAL	SAMPLE VOLUME (std.)	SAMPLE VOLUME (alt.)
cCRP	Li-Hep Plasma / Serum	5 µl	X
canine Progesterone	Li-Hep Plasma / Serum	40 µl	5 µl
bovine Progesterone	Li-Hep Plasma / Serum	40 µl	X
equine Progesterone	Li-Hep Plasma / Serum	40 µl	X
Cortisol	Li-Hep Plasma / Serum	40 µl	X
Fibrinogen	Li-Hep Plasma / Citrate Plasma	5 µl	X
Fructosamine	Li-Hep Plasma / Serum	20 µl	X
GLDH	Li-Hep Plasma / Serum	100 µl	X
Lactate	Li-Hep Plasma	20 µl	X
Lipase (pancreatic)	Li-Hep Plasma / Serum	20 µl	5 µl
NH3 Ammonia	Li-Hep Plasma	40 µl	X
SAA	Li-Hep Plasma / Serum	5 µl	X
T4	Li-Hep Plasma / Serum	20 µl	5 µl