

VET test information sheet

PARAMETER	SAMPLE MATERIAL	SAMPLE VOLUME (std.)	SAMPLE VOLUME (alt.)
cCRP	Li-hep plasma / serum / EDTA plasma	5 µl	X
canine Progesterone	Serum / li-hep plasma	40 µl	X
Fibrinogen	Li-hep plasma / citrated plasma	5 µl	X
Fructosamine	Li-hep plasma / serum / EDTA plasma	20 µl	X
GLDH	Li-hep plasma / serum	100 µl	X
Haemoglobin	Whole Blood	20 µl	X
Lactate	Plasma	20 µl	X
Lipase (pancreatic)	Li-hep plasma / serum	20 µl	5 µl
NH3	EDTA plasma	20 µl	X
Phenobarbital	Li-hep plasma / serum	20 µl	X
SAA	Serum	5 µl	X
T4	Li-hep plasma / serum / citrated plasma	20 µl	5 µl
Total Bile Acids	Serum / li-hep plasma	100 µl	X